

## To rent a CPM Unit:

- (( Call Remington Medical CPM Rental Division at **905-470-7790** or **1-800-267-5822 ext. 234** or email **cpm@remingtonmedical.com**
- (( Must provide all information including the unit required, length of time rental required, as well as any WSIB information or payment method. *Please note: all information will be kept confidential.*
- (( Remington Medical will arrange delivery and setup. Please contact us at least 3 days in advance to ensure unit availability.
- (( At the end of your rental, contact Annina in the CPM Rental Division to make arrangements for the return of the unit and to discontinue rental charges.



**Remington Medical**

*Providing excellence in home healthcare & rehabilitation products for over 30 years!*



**Remington Medical**

#9 - 401 Bentley Street, Markham, ON L3R 9T2  
t: 905-470-7790 ff: 1-800-267-5822 f: 905-470-7787 ff: 1-866-470-7787  
cpm@remingtonmedical.com www.cpmrental.ca



**Remington Medical**

# CPM

Continuous Passive Motion  
RENTAL PROGRAM

**Your First Choice for  
CPM rental needs  
in Canada.**



[www.cpmrental.ca](http://www.cpmrental.ca)

## What is CPM?

**Continuous Passive Motion (CPM)** is a postoperative, therapeutic modality that passively - without patient effort - moves a joint through a prescribed range of motion (ROM) for an extended period of time. There are several CPM devices specifically for the hand, wrist, forearm, elbow, shoulder, knee and the ankle.

This continuous motion results in significant benefits for the patient by:

- )) Helping to prevent joint stiffness and increase joint ROM
- )) Promoting the healing of the joint surfaces and soft tissue
- )) Reducing edema and pain, resulting in improved compliance with physical therapy and exercise

These physiological benefits result in a faster recovery and return to function, as well as lowering overall rehabilitation costs.

## Why Chose Remington?

- )) If you are currently using a CPM device in the hospital, we would be pleased to make arrangements for an available CPM unit for home use after discharge.
- )) We provide delivery and setup by trained professionals, at the hospital or at your home. *(Charges may apply)*
- )) We have therapists on staff to answer all your questions and concerns.
- )) We will obtain all required authorization for insurance and WSIB claims.
- )) Our rental units are thoroughly cleaned and disinfected before shipping with new soft goods.
- )) We believe in timely customer service and issues are dealt promptly.



## Our CPM Units

### E3 (Elbow)

Designed to offer the maximum in patient comfort and superior alignment accuracy. It operates at a specific speed between the programmed range of motion (ROM) extremes and is ideal for immediate post-operative use.

### WaveFlex (Hand)

Achieves full composite first flexion, plus full extension and with sophisticated finger attachments, assist each finger through a natural path of motion

### L4/L4D (Lower Limb)

Achieves a wide range of motion from -10° extension to 135° flex the and provides simultaneous hip, knee and adjustable ankle motion for total limb CPM.

### A3 (Ankle)

Provides anatomically correct motion to both the ankle and subtalar joints. Fully adjustable range of motion facilitates inversion/eversion, plantar/dorsiflexion, or a combination of both.

### S3 (Shoulder)

Designed to meet the complex demands of a shoulder CPM by providing several treatment options to achieve sequential motion in internal/external rotation and flexion/extension or abduction/adduction

### PS1 (Forearm)

delivers proven benefits to help patients regain post-operative range of motion by gently moving the forearm through a controlled 90° pronation/supination range of motion.

### W2 (Wrist)

Designed to deliver maximum mobility to the wrist joint in flexion/extension and ulnar/radial deviation or a combination of both.

### T1 (Toe)

Offers patients the proven benefits of CPM to the metatarsal head of the toe, providing easy ROM adjustments from 0° to 60° of dorsiflexion and 0° to 45° of plantar flexion.

Some Units come complete with drive unit, mounting frame, power supply, motion controller, battery charger, battery pack, carrying case and instruction manuals.

